

| MIL/MOA | MIL reticle | | | MOA reticle | | |
|---------|--|--|--|--|---|---|
| Reticle | FML-TR1 | FML-3 | FML-MT | FMA-TR1 | FMA-3 | FMA-MT |
| 4-40x52 | Line: 0.03MIL (after 2MIL: 0.06MIL) Center dot: 0.075MIL | Line: 0.03MIL (after 2MIL: 0.06MIL) Center dot: 0.075MIL | x | x | Line: 0.125MOA (about 0.04MIL) Center dot: 0.25MOA (about 0.07MIL) | x |
| 6-60x56 | Line: 0.025MIL (after 2MIL: 0.05MIL) Center dot: 0.06MIL | Line: 0.025MIL (after 2MIL: 0.05MIL) Center dot: 0.06MIL | Line: 0.02MIL Illumi center dot: 0.1MIL Non-illumi center dot: 0.035MIL | Line: 0.125MOA (about 0.04MIL) Center dot: 0.2MOA (about 0.06MIL) | Line: 0.125MOA (about 0.04MIL) Center dot: 0.2MOA (about 0.06MIL) | Line: 0.07MOA (about 0.02MIL) Center dot: 0.2MOA (about 0.06MIL) |

*The subtension of the reticle is the same either it is illuminated or non-illuminated except for the FML-MT reticle.